



## Nicolet Pocket-Dop II

- Obstetric and vascular
- ✓ Built-in Speaker
- Rechargeable Batteries



Nicolet

# Nicolet Pocket-Dop II

## Obstetric and Vascular Versatility

Pocket-Dop II is a full-featured Doppler with interchangeable obstetric and vascular probes. Use the 2 or 3 MHz OB probes to monitor the fetal heartbeat. Use the 5 or 8 MHz vascular probes to screen patients at risk for stroke, peripheral arterial disease and other occlusive vascular conditions or to monitor difficult pulses and obtain systolic blood pressures quickly and easily.

### Obstetric Examination

The Pocket-Dop II is a valuable aid in detecting pregnancy and monitoring the fetal heartbeat. The built-in speaker amplifies fetal sounds to reassure the expectant mother. A complete obstetric system should include a 3 MHz probe for early fetal heartbeat detection and a 2 MHz probe for use in later pregnancy or with larger patients. A 2 MHz waterproof probe is ideal for underwater labor or delivery.

### Vascular Applications

Use the Pocket-Dop II to compare simple Doppler blood pressures in the arm and ankle to help diagnose arterial disease in the lower extremities. Or quickly screen for lower extremity vascular conditions such as valvular incompetence.

### Emergencies and Challenging Clinical Situations

The Pocket-Dop II can be used to accurately monitor pulses and systolic blood pressures, even in the noisiest environments; confirm blood flow to wound sites, locate vessels for catheterization, withdrawal, injection or IV therapy; clearly hear low systolic blood pressures with patients who are diabetic, elderly or in shock.

### Catalog Numbers

Pocket-Dop II with:	
2 MHz probe .....	P200
3 MHz probe .....	P300
5 MHz probe .....	P500
8 MHz probe .....	P800
2 & 3 MHz probes .....	P230
2 & 5 MHz probes .....	P250
2 & 8 MHz probes .....	P280
3 & 5 MHz probes .....	P350
3 & 8 MHz probes .....	P380
5 & 8 MHz probes .....	P580
2 MHz waterproof probe .....	P2WP

### Technical Specifications

Weight:	0.285 kg (10 oz.)
Dimensions:	3.17 x 5.71 x 10.79 cm (1.25 x 2.25 x 4.25 in.)
Power Source:	110V or 220V recharger 2 replaceable AA NIMH batteries
Warranty:	5 years parts 1 year labor

This product is not CE marked



**Distributed BY:**

**GLOBAL MEDICAL DISTRIBUTION LTD**

**Website: [www.gmdltd.com](http://www.gmdltd.com) • e-mail: [info@gmdltd.com](mailto:info@gmdltd.com)**