

PODO MAX

Back to Basics

continues with the *Pod***MAX**-
a dedicated podiatry processor
to MAXimize your productivity



Over 140,000
processors sold
worldwide

*The highest throughput
dedicated table-top
podiatry processor
available anywhere!*



HOPE

PodomaX SPECIFICATIONS

The Back to Basics features of the **PodomaX** are designed to deliver high-output for its compact size and is value priced to ideally serve the moderate to high film developing requirements of the podiatry market.

Compare for yourself . . . then come to HOPE.

PodomaX versatile features:

- PodomaX utilizes Hope's unique low torque, no crossover hard roller transport system.
- PodomaX's quiet, powerful, low energy dryer easily dries film under any condition.
- PodomaX's "dual" replenishment system: may replenish either automatically or manually by push button back-up.
- Automatic energy and water saving standby system with built-in jog feature. Manual operation back-up!
- PodomaX may be set up at different speeds from 90 to 127 seconds
- Plugs into any standard wall outlet
- Cold water operation - no hot water required!
- Separate drains provide for quick draining of solutions, closed loop operation and easy use of silver recovery systems
- Light tight feed cover reduces darkroom time and eliminates costly darkroom passage systems
- Digital display of developer and dryer temperatures
- Sales and service supported by factory trained world-wide dealer network

Options available include...

- Closed loop system (no external plumbing needed)
- 8 Gallon replenishment tanks available (30.4 ltr)
- Base stand

Affordably priced - the PodomaX provides more features and benefits for less!

PodomaX

Processing Time: (Leading edge to leading edge)*

Standard Processing Cycle - 127 Seconds

Rapid Processing Cycle - 90 Seconds

* Other processing cycles available

Processing Capacity:

Standard Process:

150 (8" x 10") films/hour; 100 (10" x 12") films/hour

Rapid Process:

210 (8" x 10") films/hour; 140 (10" x 12") films/hour

Processor Tank Capacities:

Developer 1 gal/3.8 ltr

Fixer 1 gal/3.8 ltr

Wash 1 gal/3.8 ltr

Electrical:

110V, 60 Hz, 9 amp

220V, 50/60 Hz, 4.5 amp - 1 ph

Water Requirements: (GPM)

Operating 0.25 GPM/.95 ltr

Standby 0.00 gal/0.0 ltr

Dimensions:

Length (overall) 32" (81.3cm)

Width 18" (45.7cm)

Height (without base) 20" (50.8cm)

Height (with base) 46" (116.8cm)

Shipping Weight: 100 lb/45 kg

Standard:

Includes light tight feed cover, manual, and water supply hose.

Warranty:

